

ARG21851 Mouse anti-Human IgE (Fc) antibody [B3102E8] (FITC)

Package: 250 µg
Store at: 4°C

Summary

| | |
|---------------------|---|
| Product Description | FITC-conjugated Mouse Monoclonal antibody [B3102E8] recognizes Human IgE (Fc) |
| Tested Reactivity | Hu |
| Tested Application | ELISA, FACS, FLISA, ICC/IF, IHC-Fr, WB |
| Host | Mouse |
| Clonality | Monoclonal |
| Clone | B3102E8 |
| Isotype | IgG1, kappa |
| Target Name | IgE (Fc) |
| Species | Human |
| Conjugation | FITC |

Application Instructions

| Application table | Application | Dilution |
|-------------------|-------------|-----------------|
| | ELISA | Assay-dependent |
| | FACS | Assay-dependent |
| | FLISA | 1:100 - 1:400 |
| | ICC/IF | Assay-dependent |
| | IHC-Fr | Assay-dependent |
| | WB | Assay-dependent |

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

| | |
|---------------------|--|
| Form | Liquid |
| Buffer | PBS and 0.1% Sodium azide. |
| Preservative | 0.1% Sodium azide |
| Concentration | 0.5 mg/ml |
| Storage instruction | Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |