

## ARG21854 Mouse anti-Human IgE (Fc) antibody [B3102E8] (Biotin)

Package: 250 µg  
Store at: 4°C

### Summary

Product Description	Biotin-conjugated Mouse Monoclonal antibody [B3102E8] recognizes Human IgE (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, FACS, FLISA, ICC/IF, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	B3102E8
Isotype	IgG1, kappa
Target Name	IgE (Fc)
Antigen Species	Human
Conjugation	Biotin

### Application Instructions

Application table	Application	Dilution
	ELISA	1:5000 - 1:10000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.