

## ARG21876 Mouse anti-Human IgE (Fc) antibody [HP6061]

Package: 250 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [HP6061] recognizes Human IgE (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, Puri
Host	Mouse
Clonality	Monoclonal
Clone	HP6061
Isotype	IgM, kappa
Target Name	IgE (Fc)
Antigen Species	Human
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	2.5 - 10 µg/ml
	ELISPOT	Assay-dependent
	Puri	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.