

ARG21878 Mouse anti-Human IgE (Fc) antibody [HP6029] (FITC)

Package: 250 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [HP6029] recognizes Human IgE (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FLISA, Puri, WB
Host	Mouse
Clonality	Monoclonal
Clone	HP6029
Isotype	IgG2a, kappa
Target Name	IgE (Fc)
Species	Human
Conjugation	FITC

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ELISPOT	Assay-dependent
	FLISA	1:100 - 1:400
	Puri	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.