

ARG21879 Mouse anti-Human IgE (Fc) antibody [HP6029] (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Mouse Monoclonal antibody [HP6029] recognizes Human IgE (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FLISA, Puri, WB
Host	Mouse
Clonality	Monoclonal
Clone	HP6029
Isotype	IgG2a, kappa
Target Name	IgE (Fc)
Species	Human
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	ELISA	1:100 - 1:400
	ELISPOT	Assay-dependent
	FLISA	Assay-dependent
	Puri	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.