

ARG21949 Mouse anti-Human IgG1 (Fc) antibody [HP6001]

Package: 250 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [HP6001] recognizes Human IgG1 (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, FACS, FLISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	HP6001
Isotype	IgG2b, kappa
Target Name	IgG1 (Fc)
Antigen Species	Human
Immunogen	Human IgG1 myeloma protein
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	< 1 ug/ml
	FACS	Assay-dependent
	FLISA	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.