

ARG21957 Goat anti-Rat IgG antibody, pre-adsorbed

Package: 1 mg
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes Rat IgG
Tested Reactivity	Rat
Tested Application	ELISA, FACS, FLISA, ICC/IF, IHC-Fr, WB, Ifa
Specificity	The antibody reacts with the heavy chain of Rat IgG. The antibody is pre-adsorbed with Rat IgM, so the antibody may not react with Rat IgM, but may react with immunoglobulins from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG
Species	Rat
Conjugation	Un-conjugated

Application Instructions

Pre Adsorbed	Rat IgM.																	
Application table	<table><thead><tr><th>Application</th><th>Dilution</th></tr></thead><tbody><tr><td>ELISA</td><td>Assay-dependent</td></tr><tr><td>FACS</td><td>Assay-dependent</td></tr><tr><td>FLISA</td><td>Assay-dependent</td></tr><tr><td>ICC/IF</td><td>Assay-dependent</td></tr><tr><td>IHC-Fr</td><td>Assay-dependent</td></tr><tr><td>WB</td><td>Assay-dependent</td></tr><tr><td>Ifa</td><td>Assay-dependent</td></tr></tbody></table>		Application	Dilution	ELISA	Assay-dependent	FACS	Assay-dependent	FLISA	Assay-dependent	ICC/IF	Assay-dependent	IHC-Fr	Assay-dependent	WB	Assay-dependent	Ifa	Assay-dependent
Application	Dilution																	
ELISA	Assay-dependent																	
FACS	Assay-dependent																	
FLISA	Assay-dependent																	
ICC/IF	Assay-dependent																	
IHC-Fr	Assay-dependent																	
WB	Assay-dependent																	
Ifa	Assay-dependent																	
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.																	

Properties


Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	BBS (pH 8.2)
Concentration	2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw

cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.

Images



ARG21957 Goat anti-Rat IgG antibody (pre-adsorbed) ELISA image

ELISA: The plate was coated with purified Rat IgG and IgM. Immunoglobulins were detected with serially diluted ARG21957 Goat anti-Rat IgG antibody (pre-adsorbed) followed by [ARG23857](#) Swine anti-Goat IgG (H+L) antibody (HRP) (pre-adsorbed).