

## ARG22037 anti-CD45RB antibody [C363.16A] (PE-Cyanine 5)

Package: 50 µg  
Store at: 4°C

### Summary

|                     |  |
|---------------------|--|
| Product Description | PE-Cyanine 5-conjugated Rat Monoclonal antibody [C363.16A] recognizes CD45RB   |
| Tested Reactivity   | Ms   |
| Tested Application  | FACS   |
| Specificity         | Mouse CD45RB   |
| Host                | Rat  |
| Clonality           | Monoclonal   |
| Clone               | C363.16A   |
| Isotype             | IgG2a, kappa   |
| Target Name         | CD45RB   |
| Species             | Mouse  |
| Immunogen           | Cloned Mouse Th cell line  |
| Conjugation         | PE-Cyanine 5   |
| Alternate Names     | LY5; GP180; Receptor-type tyrosine-protein phosphatase C; CD45; L-CA; CD antigen CD45; Leukocyte common antigen; CD45R; LCA; T200; EC 3.1.3.48; B220 |

### Application Instructions

| Application table | <table> <tr> <th>Application</th><th>Dilution</th></tr> <tr> <td>FACS</td><td>&lt; 0.2 µg/10<sup>6</sup> cells</td></tr> </table>          | Application | Dilution | FACS | < 0.2 µg/10 <sup>6</sup> cells |
|-------------------|--|-------------|----------|------|--------------------------------|
| Application       | Dilution   |             |          |      |                                |
| FACS              | < 0.2 µg/10 <sup>6</sup> cells   |             |          |      |                                |
| Application Note  | * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist. |             |          |      |                                |


### Properties

|                     |  |
|---------------------|--|
| Form                | Liquid   |
| Buffer              | PBS, 0.1% Sodium azide and Sucrose.  |
| Preservative        | 0.1% Sodium azide  |
| Stabilizer          | Sucrose  |
| Concentration       | 0.1 mg/ml  |
| Storage instruction | Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.   |

## Bioinformation

|                |  |
|----------------|--|
| Database links | <a href="#">GeneID: 19264 Mouse</a>  |
| Gene Symbol    | PTPRC  |
| Gene Full Name | protein tyrosine phosphatase, receptor type, C   |
| Background     | <p>CD45 is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitosis, and oncogenic transformation. This PTP contains an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus is classified as a receptor type PTP. This PTP has been shown to be an essential regulator of T- and B-cell antigen receptor signaling. It functions through either direct interaction with components of the antigen receptor complexes, or by activating various Src family kinases required for the antigen receptor signaling. This PTP also suppresses JAK kinases, and thus functions as a regulator of cytokine receptor signaling. Alternatively spliced transcripts variants of this gene, which encode distinct isoforms, have been reported. [provided by RefSeq, Jun 2012]</p> |
| Function       | <p>CD45: Protein tyrosine-protein phosphatase required for T-cell activation through the antigen receptor. Acts as a positive regulator of T-cell coactivation upon binding to DPP4. The first PTPase domain has enzymatic activity, while the second one seems to affect the substrate specificity of the first one. Upon T-cell activation, recruits and dephosphorylates SKAP1 and FYN. Dephosphorylates LYN, and thereby modulates LYN activity.</p> <p>(Microbial infection) Acts as a receptor for human cytomegalovirus protein UL11 and mediates binding of UL11 to T-cells, leading to reduced induction of tyrosine phosphorylation of multiple signaling proteins upon T-cell receptor stimulation and impaired T-cell proliferation. [UniProt]</p>   |
| Research Area  | Developmental Biology antibody; Immune System antibody; Neuroscience antibody; Signaling Transduction antibody; Mouse Inflammatory Cell Marker antibody; B Cell Marker antibody  |
| Calculated Mw  | 147 kDa  |
| PTM            | Heavily N- and O-glycosylated.   |

## Images



ARG22037 anti-CD45RB antibody [C363.16A] (PE-Cyanine 5) FACS image

Flow Cytometry: BALB/c Mouse splenocytes stained with ARG22037 anti-CD45RB antibody [C363.16A] (PE-Cyanine 5).