

ARG22117 anti-MHC Class I H2 Db antibody [27-11-13S]

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [27-11-13S] recognizes MHC Class I H2 Db
Tested Reactivity	Ms
Tested Application	FACS, ICC/IF, IP
Specificity	Mouse H-2Db
Host	Mouse
Clonality	Monoclonal
Clone	27-11-13S
Isotype	IgG2a, kappa
Target Name	MHC Class I H2 Db
Species	Mouse
Immunogen	BDF1 Mouse splenocytes
Conjugation	Un-conjugated
Alternate Names	H-2 class I histocompatibility antigen, D-B alpha chain; H2-DB; H-2D

Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
	ICC/IF	Assay-dependent
	IP	Assay-dependent
	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Application Note		

Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Database links	GeneID: 14964 Mouse Swiss-port # P01899 Mouse
Gene Symbol	H2-D1
Gene Full Name	histocompatibility 2, D region locus 1
Function	Involved in the presentation of foreign antigens to the immune system. [UniProt]
Calculated Mw	41 kDa