

ARG22684 anti-Candida albicans antibody (Biotin)

Package: 500 µl
Store at: 4°C

Summary

Product Description	Biotin-conjugated Rabbit Polyclonal antibody recognizes Candida albicans This antibody recognizes numerous proteins in a soluble C. albicans extract in immunoelectrophoresis. Candida albicans is the most frequently isolated fungal pathogen of Humans. It principally infects the skin and mucosal surfaces causing candidiasis or thrush. This antibody has not been absorbed and cross reacts with other yeasts.
Tested Reactivity	C. albicans
Tested Application	ELISA, ICC/IF, IHC-Fr
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Candida albicans
Species	Yeast
Immunogen	Candida albicans, type A
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ICC/IF	1:200
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Unpurified.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	4 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.