

ARG23773 Goat IgG Isotype Control antibody (Spectral Red)

Package: 100 tests
Store at: 4°C

Summary

Product Description	Spectral Red-conjugated Goat Polyclonal antibody as a negative control antibody for Goat IgG
Tested Application	BL, ELISA, FACS, FLISA, IHC-Fr, IHC-P
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Goat IgG
Conjugation	Spectral Red

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	Assay-dependent
	FACS	10 µl/10 ⁶ cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.