

ARG23844 Mouse anti-Rabbit Light Chain antibody [SB62a] (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Mouse Monoclonal antibody [SB62a] recognizes Rabbit Light Chain
Tested Reactivity	Rb
Tested Application	ELISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	SB62a
Isotype	IgG2b
Target Name	Light Chain
Species	Rabbit
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	ELISA	1:2000 - 1:4000
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol.
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.