

ARG23932 Mouse IgG1 Isotype Control antibody [15H6] (FITC)

Package: 100 tests
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [15H6] as a negative control antibody for Mouse IgG1
Tested Application	BL, FACS, FLISA, ICC/IF, IHC-Fr, IHC-P, WB
Specificity	T-2 mycotoxin
Host	Mouse
Clonality	Monoclonal
Clone	15H6
Isotype	IgG1, kappa
Target Name	Mouse IgG1
Conjugation	FITC


Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	FACS	10 µl/10 ⁶ cells
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG23932 Mouse IgG1 Isotype Control antibody [15H6] (FITC) FLISA image

FLISA: The plate was coated with [ARG21883](#) Goat anti-Mouse IgG1 antibody (pre-adsorbed), [ARG23801](#) Goat anti-Mouse IgG2a antibody (pre-adsorbed), [ARG21891](#) Goat anti-Mouse IgG2b antibody (pre-adsorbed), Goat anti-Mouse IgG3 antibody, [ARG21539](#) Goat anti-Mouse IgA antibody (pre-adsorbed) and [ARG21517](#) Goat anti-Mouse IgM antibody (pre-adsorbed). Serially diluted ARG23932 Mouse IgG1 Isotype Control antibody [15H6] (FITC) was captured and fluorescence intensity quantified.