

ARG30281 PLA2G7 / Lp-PLA2 ELISA Antibody Duo

Package: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10615	anti-PLA2G7 / Lp-PLA2 antibody [PL42cc]	Mouse mAb	Hu	ELISA, WB	50 µg
ARG10616	anti-PLA2G7 / Lp-PLA2 antibody [PL46cc]	Mouse mAb	Hu	ELISA, WB	50 µg

Summary

Target Name	PLA2G7 / Lp-PLA2
Alternate Names	PLA2G7 / Lp-PLA2 ELISA antibody; PLA2G7 / Lp-PLA2 antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Full Name	ELISA Antibody Duo for PLA2G7 / Lp-PLA2
Highlight	Related Product: anti-PLA2G7 / Lp-PLA2 antibody;