

ARG30289 TSH / Thyroid Stimulating Hormone ELISA Antibody Duo

Package: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10635	anti-TSH / Thyroid Stimulating Hormone antibody [1CT1]	Mouse mAb	Hu	ELISA, WB, sELISA	50 µg
ARG10636	anti-TSH / Thyroid Stimulating Hormone antibody [11E4cc]	Mouse mAb	Hu	ELISA, sELISA	50 µg

Summary


Target Name	TSH / Thyroid Stimulating Hormone
Alternate Names	TSH / Thyroid Stimulating Hormone ELISA antibody; Thyroid Stimulating Hormone ELISA antibody; TSH ELISA antibody; TSH / Thyroid Stimulating Hormone antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
---------------------	--


Note For laboratory research only, not for drug, diagnostic or other use.

Images




ARG10635 anti-TSH / Thyroid Stimulating Hormone antibody [1CT1]
ELISA image

ELISA: Titration curves of ARG10635 anti-TSH / Thyroid Stimulating Hormone antibody [1CT1]. Antigens: Canine recombinant TSH (in-house preparation) and Human native TSH (NHPP), 0.03 µg/well.



ARG10635 anti-TSH / Thyroid Stimulating Hormone antibody [1CT1] sELISA image

Sandwich ELISA: Calibration curves for MAb combinations 11E4cc - 1CT1. The immunoassay was performed in 96-well streptavidin-coated plates (Greiner Bio-One). Capture antibodies: ARG10636 anti-Thyroid Stimulating Hormone (TSH) antibody [11E4cc] (additional biotinylated preparation). Detection antibody: ARG10635 anti-TSH / Thyroid Stimulating Hormone antibody [1CT1] (Eu3+-labeled). Antigens: Canine recombinant TSH (in-house preparation) and native Human TSH (NHPP).



ARG10636 anti-TSH / Thyroid Stimulating Hormone antibody [11E4cc] ELISA image

ELISA: Titration curves of ARG10636 anti-TSH / Thyroid Stimulating Hormone antibody [11E4cc]. Antigens: Canine recombinant TSH (in-house preparation) and Human native TSH (NHPP), 0.03 µg/well.