

## ARG30336 Human MDSC Marker Antibody Duo

Package: 1 pair  
Store at: -20°C

### Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG22000	anti-CD11b antibody [M1/70]	Rat mAb	Hu, Ms, Bb, R. Mk	BL, FACS, ICC/IF, IHC-Fr, IP	50 µg
ARG65560	anti-CD33 antibody [WM53]	Mouse mAb	Hu, NHuPrm	CytoF®-candidate, FACS, FuncSt, ICC/IF, IHC-Fr, IP, WB	50 µg

### Summary

**Product Description** Human MDSC Marker Antibody Duo includes antibodies to MDSC marker, CD11b and CD33, for Myeloid-derived suppressor cells (MDSCs) study. Myeloid-derived suppressor cells (MDSCs) are a heterogeneous population of immune cells from the myeloid lineage. MDSCs are well-known key negative regulators of the immune response in tumor microenvironment. Accumulated evidences show that MDSCs may serve as a therapeutic target for preventing tumor progression. In human, MDSCs are defined as cells expressing CD33 and CD11b. arigo's Human MDSC Marker Antibody Duo comprises CD33 and CD11b antibodies. It is the best solution to identify the MDSCs in human tumor tissues.

Related news:

[New antibody panels and duos for Tumor immune microenvironment](#)  
[Anti-SerpinB9 therapy, a new strategy for cancer therapy](#)

**Target Name** Human MDSC marker

**Alternate Names** Human MDSC Marker antibody; Myeloid-derived suppressor cells marker antibody; MDSCs marker antibody; CD11b antibody; CD33 antibody

### Properties

**Storage instruction** For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.


**Note** For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

**Gene Full Name** Human Myeloid-derived suppressor cell (MDSC) marker Antibody Duo


**Highlight** Related products:  
[anti-CD11b antibody](#);  
[anti-CD33 antibody](#);  
[MDSC Duos / Panels](#);

**Research Area** Myeloid-derived suppressor cells marker; MDSCs marker; Tumor immune microenvironment; tumor-associated myeloid cells




ARG22000 anti-CD11b antibody [M1/70] IHC-Fr image

Immunohistochemistry: Cryosection of Mouse lymph node stained with ARG22000 anti-CD11b antibody [M1/70] (green, A) and anti-Mouse CD8 antibody (red, B). Merged image in C.




ARG65560 anti-CD33 antibody [WM53] IHC-Fr image

Immunohistochemistry: Frozen section of Human colon tissue stained with ARG65560 anti-CD33 antibody [WM53].




ARG65560 anti-CD33 antibody [WM53] CyTOF image

CyTOF: Human peripheral blood stained with ARG65560 anti-CD33 antibody [WM53] (149Sm). Singlet cells were gated for data analysis.



ARG65560 anti-CD33 antibody [WM53] FACS image

Flow Cytometry: Separation of Human CD33 positive Monocytes (red) from Human CD33 negative Lymphocytes (black-dashed). Human peripheral blood stained with ARG65560 anti-CD33 antibody [WM53], followed by incubation with APC labelled Goat anti-Mouse secondary antibody.



#### ARG65560 anti-CD33 antibody [WM53] FACS image

Flow Cytometry: PBMC stained with ARG65560 anti-CD33 antibody [WM53] at  $0.5 \mu\text{g}/10^6$  cells (right histogram) or isotype control (left histogram), followed by incubation with PE labelled secondary antibody.