

ARG30338 Multiplex Secondary Antibody Panel

Package: 1 kit
Store at: -20°C, 4°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG65350	Goat anti-Mouse IgG antibody (HRP)	Goat pAb	Ms	ELISA, IHC-P, WB	20 µl
ARG65351	Goat anti-Rabbit IgG antibody (HRP)	Goat pAb	Rb	ELISA, IHC-P, WB	20 µl
ARG65352	Donkey anti-Goat IgG antibody (HRP)	Donkey pAb	Goat	ELISA, IHC-P, WB	20 µl
ARG21747	Goat anti-Chicken IgY (H+L) antibody (HRP)	Goat pAb	Chk	ELISA, FACS, FLISA, ICC/IF, IHC-P, WB	20 µl
ARG21771	Rabbit anti-Rat IgG (H+L) antibody (HRP)	Rabbit pAb	Rat	ELISA, FLISA, IHC-Fr, IHC-P, WB	20 µl

Summary


Target Name Multiplex secondary antibody

Properties

Storage instruction Store antibodies at 4°C or -20°C. Please refer to the each product datasheet for detail temperatures of the antibodies.

Note For laboratory research only, not for drug, diagnostic or other use.


Images



ARG65350 Goat anti-Mouse IgG antibody (HRP) WB image

Western blot: Rat basolateral amygdala stained with [ARG62347 anti-beta Tubulin antibody \[BT7R\]](#) at 1:1000 dilution, ARG65350 Goat anti-Mouse IgG antibody (HRP) at 1:5000 dilution.

From Guang-Bing Duan et al. CNS Neurosci Ther. (2024), [doi: 10.1111/cns.14611](https://doi.org/10.1111/cns.14611), Fig. 4.D.



ARG65350 Goat anti-Mouse IgG antibody (HRP) IHC-P image


From Cheng-Feng Chu et al. J Pers Med. (2021), [doi: 10.3390/jpm11121326](https://doi.org/10.3390/jpm11121326), Fig. 6.



ARG65351 Goat anti-Rabbit IgG antibody (HRP) WB image


Western blot: Mouse retina stained with [ARG65693 anti-alpha Tubulin antibody](#) and ARG65351 Goat anti-Rabbit IgG antibody (HRP)

From Xiaoyuan Ye et al. Mol Ther Nucleic Acids. (2024), [doi: 10.1016/j.omtn.2024.102209](https://doi.org/10.1016/j.omtn.2024.102209), Fig. 5.D.



ARG65351 Goat anti-Rabbit IgG antibody (HRP) IHC-P image

From Yu-Qian Song et al. J Mol Med (Berl) (2022), [doi: 10.1007/s00109-021-02165-0](https://doi.org/10.1007/s00109-021-02165-0), Fig. 5.c.



ARG65352 Donkey anti-Goat IgG antibody (HRP) WB image

Western blot: Endometrial Stromal cells stained with [ARG63922 anti-EP4 prostaglandin Receptor antibody](#) at 1:1000 dilution and ARG65352 Donkey anti-Goat IgG antibody (HRP).

From Qingqing Huang et al. Reprod Med Biol. (2021), [doi: 10.1002/rmb2.12423](https://doi.org/10.1002/rmb2.12423), Fig. 3.C.