

ARG30360 Adipogenic Marker Antibody Panel

Package: 1 pair

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG55241	anti-PPAR gamma antibody	Rabbit pAb	Hu, Ms, Rat	WB	20 µl
ARG66798	anti-C / EBP alpha antibody	Rabbit pAb	Hu, Ms, Rat	IHC-P, WB	20 µl
ARG42921	anti-FABP4 antibody	Rabbit pAb	Hu, Ms, Rat	ICC/IF, WB	20 µl

Summary


Product Description	Adipogenic Marker Antibody Panel is an all-in-one solution to make Adipogenic Marker research easy and economic.
Target Name	PPAR gamma; C / EBP alpha; FABP4

Properties

Storage instruction	-20°C
Note	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation


Gene Full Name	Antibody Panel for Adipogenic Marker
Highlight	Related Product: anti-PPAR gamma antibody; anti-C / EBP alpha antibody; anti-FABP4 antibody;



ARG55241 anti-PPAR gamma antibody IHC-P image


Immunohistochemistry: Mouse bone stained with ARG55241 anti-PPAR gamma antibody.

From Jingyuan Wen et al. Biomed Pharmacother., (2023) [doi: 10.1016/j.biopha.2022.114122](https://doi.org/10.1016/j.biopha.2022.114122), Fig. 5C.




ARG55241 anti-PPAR gamma antibody WB image

Western blot: 30 µg of MCF7 and MDA-MB-231 (PPAR gamma low expression cell line) stained with ARG55241 anti-PPAR gamma antibody at 1:1000 dilution.




ARG55241 anti-PPAR gamma antibody WB image

Western blot: 30 µg of Rat liver lysate stained with ARG55241 anti-PPAR gamma antibody at 1:1000 dilution.



ARG66798 anti-C/EBP alpha antibody IHC-P image

Immunohistochemistry: Formalin-fixed and paraffin-embedded Human lung cancer tissue. Antigen Retrieval: Heat mediation was performed in Sodium citrate buffer (pH 6.0). The tissue section was stained with ARG66798 anti-C/EBP alpha antibody at room temperature. The section was counterstained with haematoxylin and mounted with DPX.



ARG66798 anti-C/EBP alpha antibody WB image

Western blot: Mouse embryo, Rat ovary, NIH3T3L1, PC-12, HeLa, Panc1 and SGC7901 whole cell lysates stained with ARG66798 anti-C/EBP alpha antibody.



ARG42921 anti-FABP4 antibody WB image

Western blot: Human fetal heart lysate stained with ARG42921 anti-FABP4 antibody.

Human fetal heart