

ARG30361 ABC Transporter Antibody Panel

Package: 1 pair

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG40732	anti-MDR1 / P Glycoprotein 1 antibody	Rabbit pAb	Hu, Ms, Rat	IHC-P, WB	20 µl
ARG67310	anti-ABCC1 / MRP1 antibody [6B02]	Rabbit mAb	Hu, Ms	IHC-P, WB	20 µl
ARG66193	anti-ABCG2 antibody	Rabbit pAb	Hu, Ms	ICC/IF, IHC-P, WB	20 µl

Summary


Product Description	ABC Transporter Antibody Panel is an all-in-one solution to make ABC Transporter research easy and economic.
Target Name	MDR1 / P Glycoprotein 1; ABCC1 / MRP1 ; ABCG2

Properties

Storage instruction	-20°C
Note	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation


Gene Full Name	Antibody Panel for ABC Transporter
Highlight	Related Product: anti-MDR1 / P Glycoprotein 1 antibody; anti-ABCC1 / MRP1 antibody; anti-ABCG2 antibody;



ARG40732 anti-MDR1 / P Glycoprotein 1 antibody WB image


Western blot: SOC stained with ARG40732 anti-MDR1 / P Glycoprotein 1 antibody.

From Zhang J et al. J Ovarian Res- (2021), [doi: 10.1186/s13048-021-00873-2](https://doi.org/10.1186/s13048-021-00873-2), Fig. 1. C.




ARG40732 anti-MDR1 / P Glycoprotein 1 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human liver stained with ARG40732 anti-MDR1 / P Glycoprotein 1 antibody.




ARG40732 anti-MDR1 / P Glycoprotein 1 antibody WB image

Western blot: HeLa cell lysate stained with ARG40732 anti-MDR1 / P Glycoprotein 1 antibody.



ARG67310 anti-ABCC1 / MRP1 antibody [6B02] IHC-P image

Immunohistochemistry: Human lung cancer stained with ARG67310 anti-ABCC1 / MRP1 antibody [6B02].




ARG67310 anti-ABCC1 / MRP1 antibody [6B02] WB image

Western blot: HCT-16 stained with ARG67310 anti-ABCC1 / MRP1 antibody [6B02].




ARG67310 anti-ABCC1 / MRP1 antibody [6B02] WB image

Western blot: Mouse lung stained with ARG67310 anti-ABCC1 / MRP1 antibody [6B02].




ARG66193 anti-ABCG2 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human brain tissue stained with ARG66193 anti-ABCG2 antibody at 1:100 dilution.



ARG66193 anti-ABCG2 antibody WB image

Western blot: NIH-3T3 and 293 cell lysates stained with ARG66193 anti-ABCG2 antibody at 1:500 dilution.



ARG66193 anti-ABCG2 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human brain tissue stained with ARG66193 anti-ABCG2 antibody 1:100 dilution.