

ARG40932 anti-Kininogen 1 antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes Kininogen 1
Tested Reactivity	Ms
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Kininogen 1
Species	Mouse
Immunogen	Synthetic peptide corresponding to aa. 227-259 of Mouse Kininogen 1. (ECRGNLFMDINNKIANFSQSCTLYSGDDLVEAL)
Conjugation	Un-conjugated
Alternate Names	Williams-Fitzgerald-Flaujeac factor; Kallidin II; High molecular weight kininogen; KNG; Fitzgerald factor; Alpha-2-thiol proteinase inhibitor; BDK; BK; Kininogen-1; HMWK; Kallidin I; Ile-Ser-Bradykinin

Application Instructions

Application table	Application	Dilution
	IHC-P	1:200 - 1:1000
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	0.2% Na ₂ HPO ₄ , 0.9% NaCl, 0.05% Sodium azide and 5% BSA.
Preservative	0.05% Sodium azide
Stabilizer	5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	KNG1
Gene Full Name	kininogen 1
Background	This gene uses alternative splicing to generate two different proteins- high molecular weight kininogen (HMWK) and low molecular weight kininogen (LMWK). HMWK is essential for blood coagulation and assembly of the kallikrein-kinin system. Also, bradykinin, a peptide causing numerous physiological effects, is released from HMWK. Bradykinin also functions as an antimicrobial peptide with antibacterial and antifungal activity. In contrast to HMWK, LMWK is not involved in blood coagulation. Three transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Nov 2014]
Function	(1) Kininogens are inhibitors of thiol proteases; (2) HMW-kininogen plays an important role in blood coagulation by helping to position optimally prekallikrein and factor XI next to factor XII; (3) HMW-kininogen inhibits the thrombin- and plasmin-induced aggregation of thrombocytes; (4) the active peptide bradykinin that is released from HMW-kininogen shows a variety of physiological effects: (4A) influence in smooth muscle contraction, (4B) induction of hypotension, (4C) natriuresis and diuresis, (4D) decrease in blood glucose level, (4E) it is a mediator of inflammation and causes (4E1) increase in vascular permeability, (4E2) stimulation of nociceptors (4E3) release of other mediators of inflammation (e.g. prostaglandins), (4F) it has a cardioprotective effect (directly via bradykinin action, indirectly via endothelium-derived relaxing factor action); (5) LMW-kininogen inhibits the aggregation of thrombocytes; (6) LMW-kininogen is in contrast to HMW-kininogen not involved in blood clotting. [UniProt]
Calculated Mw	72 kDa
PTM	Bradykinin is released from kininogen by plasma kallikrein. Hydroxylation of Pro-383 occurs prior to the release of bradykinin. Phosphorylated by FAM20C in the extracellular medium. N- and O-glycosylated. O-glycosylated with core 1 or possibly core 8 glycans. [UniProt]
Cellular Localization	Secreted, extracellular space. [UniProt]

Images



ARG40932 anti-Kininogen 1 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Mouse kidney tissue stained with ARG40932 anti-Kininogen 1 antibody.



ARG40932 anti-Kininogen 1 antibody WB image

Western blot: 50 µg of Mouse lung, 50 µg of Mouse testis, 50 µg of Mouse liver, 40 µg of HEPA and 40 µg of NEURO whole cell lysates stained with ARG40932 anti-Kininogen 1 antibody at 0.5 µg/ml dilution.