

# ARG43017 anti-PTPN1 / PTP1B phospho (Ser378) antibody

Package: 100 μl Store at: -20°C

## Summary

Product Description	Rabbit Polyclonal antibody recognizes PTPN1 / PTP1B phospho (Ser378)
Tested Reactivity	Hu
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	lgG
Target Name	PTPN1 / PTP1B
Species	Human
Immunogen	Phosphospecific peptide around Ser378 of Human PTPN1 / PTP1B.
Conjugation	Un-conjugated
Alternate Names	PTP1B; EC 3.1.3.48; Tyrosine-protein phosphatase non-receptor type 1; PTP-1B; Protein-tyrosine phosphatase 1B

### **Application Instructions**

Application table	Application	Dilution
	IHC-P	1:50
	WB	1:2000 - 1:10000
Application Note	* The dilutions indicate recomme should be determined by the scie	nded starting dilutions and the optimal dilutions or concentrations ntist.
Positive Control	HL-60	
Observed Size	~ 52 kDa	

### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	50 mM Tris-Glycine (pH 7.4), 150 mM NaCl, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Preservative	0.01% Sodium azide
Stabilizer	40% Glycerol and 0.05% BSA
Concentration	Batch dependent
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

#### Bioinformation

Gene Symbol	PTPN1
Gene Full Name	protein tyrosine phosphatase, non-receptor type 1
Background	The protein encoded by this gene is the founding member of the protein tyrosine phosphatase (PTP) family, which was isolated and identified based on its enzymatic activity and amino acid sequence. PTPs catalyze the hydrolysis of the phosphate monoesters specifically on tyrosine residues. Members of the PTP family share a highly conserved catalytic motif, which is essential for the catalytic activity. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP has been shown to act as a negative regulator of insulin signaling by dephosphorylating the phosphotryosine residues of insulin receptor kinase. This PTP was also reported to dephosphorylate epidermal growth factor receptor kinase, as well as JAK2 and TYK2 kinases, which implicated the role of this PTP in cell growth control, and cell response to interferon stimulation. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2013]
Function	Tyrosine-protein phosphatase which acts as a regulator of endoplasmic reticulum unfolded protein response. Mediates dephosphorylation of EIF2AK3/PERK; inactivating the protein kinase activity of EIF2AK3/PERK. May play an important role in CKII- and p60c-src-induced signal transduction cascades. May regulate the EFNA5-EPHA3 signaling pathway which modulates cell reorganization and cell-cell repulsion. May also regulate the hepatocyte growth factor receptor signaling pathway through dephosphorylation of MET. [UniProt]
Calculated Mw	50 kDa
РТМ	Oxidized on Cys-215; the Cys-SOH formed in response to redox signaling reacts with the alpha-amido of the following residue to form a sulfenamide cross-link, triggering a conformational change that inhibits substrate binding and activity. The active site can be restored by reduction.
	Ser-50 is the major site of phosphorylation as compared to Ser-242 and Ser-243. Activated by phosphorylation at Ser-50.
	S-nitrosylation of Cys-215 inactivates the enzyme activity.
	Sulfhydration at Cys-215 following endoplasmic reticulum stress inactivates the enzyme activity, promoting EIF2AK3/PERK activity. [UniProt]
Cellular Localization	Endoplasmic reticulum membrane; Peripheral membrane protein; Cytoplasmic side. Note=Interacts with EPHA3 at the cell membrane. [UniProt]

#### Images



# ARG43017 anti-PTPN1 / PTP1B phospho (Ser378) antibody WB image

Western blot: HL-60 cell lysate stained with ARG43017 anti-PTPN1 / PTP1B phospho (Ser378) antibody at 1:1000 dilution.