

## ARG45011 anti-IgG (Fc) antibody [RMG02]

Package: 50 µg  
Store at: -20°C

### Summary

Product Description	Goat monoclonal antibody [RMG02] recognizes Rabbit IgG (Fc)
Tested Reactivity	Rb
Tested Application	ELISA, WB
Host	Goat
Clonality	Monoclonal
Clone	RMG02
Isotype	IgG
Target Name	IgG (Fc)
Species	Rabbit
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	0.01 - 0.5 µg/mL
	WB	0.1 - 0.5 µg/ml
Application Note	WB: This antibody only reacts to nonreduced Rabbit IgG (Fc). * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Protein A affinity purified.
Buffer	PBS with 50% Glycerol, 1% BSA and 0.09% sodium azide
Preservative	0.09% sodium azide
Stabilizer	50% Glycerol, 1% BSA and 0.09%
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.