

ARG45034 anti-25-OH Vitamin D3 antibody antibody [RMH04]

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit monoclonal antibody [RMH04] recognizes 25-OH Vitamin D3
Tested Reactivity	Other
Tested Application	ELISA
Host	Rabbit
Clonality	Monoclonal
Clone	RMH04
Isotype	IgG
Target Name	25-OH Vitamin D3
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	50 ng/well - 200 ng/well

Application Note