

ARG46728 anti-SDHC antibody

Package: 50 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes SDHC
Tested Reactivity	Hu, Ms, Rat
Tested Application	ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	SDHC
Conjugation	Un-conjugated
Alternate Names	Succinate dehydrogenase cytochrome b560 subunit, mitochondrial; QPs1; Integral membrane protein CII-3; SDH3; CYB560; PGL3; Succinate dehydrogenase complex subunit C; Succinate-ubiquinone oxidoreductase cytochrome B large subunit; QPS1; QPs-1; CYBL

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50 - 1:200
	IHC-P	1:50 - 1:200
	WB	1:500 - 1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Observed Size	12 kDa	


Properties

Form	Liquid
Purification	Affinity chromatography purified
Buffer	PBS, Proclin 300 or sodium azide and 50% Glycerol
Preservative	Proclin 300 or sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	SDHC
Gene Full Name	succinate dehydrogenase complex, subunit C, integral membrane protein, 15kDa
Background	This gene encodes one of four nuclear-encoded subunits that comprise succinate dehydrogenase, also known as mitochondrial complex II, a key enzyme complex of the tricarboxylic acid cycle and aerobic respiratory chains of mitochondria. The encoded protein is one of two integral membrane proteins that anchor other subunits of the complex, which form the catalytic core, to the inner mitochondrial membrane. There are several related pseudogenes for this gene on different chromosomes. Mutations in this gene have been associated with paragangliomas. Alternatively spliced transcript variants have been described. [provided by RefSeq, May 2013]
Function	Membrane-anchoring subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone (coenzyme Q). [UniProt]
Calculated Mw	19 kDa
Cellular Localization	Mitochondrion inner membrane; Multi-pass membrane protein. [UniProt]

Images



ARG46728 anti-SDHC antibody IHC-P image

Immunohistochemistry: Rat kidney stained with ARG46728 anti-SDHC antibody.




ARG46728 anti-SDHC antibody ICC/IF image

Immunofluorescence: HeLa stained with ARG46728 anti-SDHC antibody.



ARG46728 anti-SDHC antibody WB image

Western blot: HeLa stained with ARG46728 anti-SDHC antibody.



ARG46728 anti-SDHC antibody WB image

Western blot: Mouse liver stained with ARG46728 anti-SDHC antibody.