

ARG46752 anti-TPD52 antibody

Package: 50 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes TPD52
Tested Reactivity	Hu, Rat
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	TPD52
Species	Human
Immunogen	Recombinant protein of Human TPD52
Conjugation	Un-conjugated
Alternate Names	PrLZ; Tumor protein D52; PC-1; Protein N8; hD52; N8L; D52

Application Instructions

Application table	Application	Dilution
	IHC-P	1:50-1:200
	WB	1:1000-1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Observed Size	26 kDa	


Properties

Form	Liquid
Purification	Affinity purification with immunogen
Buffer	PBS, 50% Glycerol, Proclin 300 or sodium azide
Preservative	Proclin 300 or sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	TPD52
Gene Full Name	tumor protein D52
Research Area	Cancer antibody; Immune System antibody
Calculated Mw	24 kDa

Images




ARG46752 anti-TPD52 antibody IHC-P image

Immunohistochemistry: Human tonsil stained with ARG46752 anti-TPD52 antibody.



ARG46752 anti-TPD52 antibody WB image

Western blot: HeLa stained with ARG46752 anti-TPD52 antibody.



ARG46752 anti-TPD52 antibody IHC-P image

Immunohistochemistry: Rat brain stained with ARG46752 anti-TPD52 antibody.