

## ARG54042 anti-DsbA tag antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody recognizes DsbA tag
Tested Reactivity	Other
Tested Application	WB
Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Target Name	DsbA tag
Species	Others
Immunogen	Purified recombinant E.coli. DsbA protein fragments over expressed in E.coli.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Observed Size	24 kDa	

### Properties


Form	Liquid
Purification	Affinity purified
Buffer	0.1M Tris-Glycine (pH 7.4), 150 mM NaCl, 0.2% Sodium azide, 0.1mg/ml BSA and 50% Glycerol
Preservative	0.2% Sodium azide
Stabilizer	0.1mg/ml BSA, 50% Glycerol
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Background	Required for disulfide bond formation in some periplasmic proteins such as PhoA or OmpA. Acts by transferring its disulfide bond to other proteins and is reduced in the process. DsbA is reoxidized by DsbB. It is required for pilus biogenesis
Research Area	Cell Biology and Cellular Response antibody; Controls and Markers antibody; Metabolism antibody; Signaling Transduction antibody
Cellular Localization	Periplasm

## Images

---



ARG54042 anti-DsbA tag antibody WB image

Western blot: DsbA Recombinant protein stained with ARG54042 anti-DsbA tag antibody at 1:1000 dilution.