

ARG54547
anti-V5 tag antibodyPackage: 50 µg
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes V5 tag
Tested Reactivity	Other
Tested Application	ELISA, ICC/IF, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	V5 tag
Species	Others
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.