

## ARG54576 anti-Secretogranin II antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes Secretogranin II
Tested Reactivity	Hu, Rat
Tested Application	IHC, IP, WB
Specificity	This antibody reacts specifically with rat and human secretogranin II.
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Secretogranin II
Species	Human
Immunogen	N-terminal sequence (amino acids 1-19) of rat/human secretogranin II coupled to ovalbumin.
Conjugation	Un-conjugated
Alternate Names	EM66; Secretogranin II; CHGC; SgII; SN; Secretogranin-2; Chromogranin-C

### Application Instructions

Application Note	<p>This antibody has been used successfully in Western blot assays (1:2,500 - 1:10,000 dilution) and immuno precipitation of secretogranin II (1:100 - 1:500 dilution), as well as in immunohistochemical identification of secretogranin II (1:1,000 dilution).</p> <p>* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.</p>
------------------	--

### Properties

Form	Liquid
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Database links	<a href="#">GeneID: 24765 Rat</a>
	<a href="#">GeneID: 7857 Human</a>
	<a href="#">Swiss-port # P10362 Rat</a>
	<a href="#">Swiss-port # P13521 Human</a>

Gene Symbol	SCG2
Gene Full Name	secretogranin II
Background	The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Studies in rodents suggest that the full-length protein, secretogranin II, is involved in the packaging or sorting of peptide hormones and neuropeptides into secretory vesicles. The full-length protein is cleaved to produce the active peptide secretoneurin, which exerts chemotaxic effects on specific cell types, and EM66, whose function is unknown. [provided by RefSeq, Jul 2008]
Function	Secretogranin-2 is a neuroendocrine secretory granule protein, which is the precursor for biologically active peptides. [UniProt]
Research Area	Neuroscience antibody
Calculated Mw	71 kDa
PTM	O-glycosylated.