

ARG58295 anti-Argininosuccinate Lyase antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes Argininosuccinate Lyase
Tested Reactivity	Ms, Rat
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Argininosuccinate Lyase
Species	Human
Immunogen	Synthetic peptide corresponding to a sequence of Human Argininosuccinate Lyase (YTHLQRAQPIRWSHWILSHAVALTDRSERLLEVRKRIN).
Conjugation	Un-conjugated
Alternate Names	EC 4.3.2.1; Argininosuccinate lyase; ASAL; Arginosuccinase

Application Instructions

Application table	Application	Dilution
	WB	0.1 - 0.5 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.05% Sodium azide and 5% BSA.
Preservative	0.05% Sodium azide
Stabilizer	5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	ASL
Gene Full Name	argininosuccinate lyase
Background	<p>This gene encodes a member of the lyase 1 family. The encoded protein forms a cytosolic homotetramer and primarily catalyzes the reversible hydrolytic cleavage of argininosuccinate into arginine and fumarate, an essential step in the liver in detoxifying ammonia via the urea cycle. Mutations in this gene result in the autosomal recessive disorder argininosuccinic aciduria, or argininosuccinic acid lyase deficiency. A nontranscribed pseudogene is also located on the long arm of chromosome 22. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]</p>
Calculated Mw	52 kDa
PTM	<p>Acetylation modifies enzyme activity in response to alterations of extracellular nutrient availability. Acetylation increased with trichostin A (TSA) or with nicotinamide (NAM). Glucose increases acetylation by about a factor of 3 with decreasing enzyme activity. Acetylation on Lys-288 is decreased on the addition of extra amino acids resulting in activation of enzyme activity. [UniProt]</p>

Images



ARG58295 anti-Argininosuccinate Lyase antibody WB image

Western blot: 50 µg of Rat liver, Rat kidney, Rat lung and Mouse liver lysates stained with ARG58295 anti-Argininosuccinate Lyase antibody at 0.5 µg/ml dilution, overnight at 4°C, under reducing conditions.