

ARG59584 anti-HAS3 antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes HAS3
Tested Reactivity	Hu, Ms, Rat
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	HAS3
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 67-281 of Human HAS3 (NP_619515.1).
Conjugation	Un-conjugated
Alternate Names	HA synthase 3; EC 2.4.1.212; Hyaluronic acid synthase 3; Hyaluronate synthase 3; Hyaluronan synthase 3

Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Rat heart, Mouse heart and SKOV3	
Observed Size	73 kDa	


Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	HAS3
Gene Full Name	hyaluronan synthase 3
Background	The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyaluronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene family. Compared to the proteins encoded by other members of this gene family, this protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010]
Function	Catalyzes the addition of GlcNAc or GlcUA monosaccharides to the nascent hyaluronan polymer. Therefore, it is essential to hyaluronan synthesis a major component of most extracellular matrices that has a structural role in tissues architectures and regulates cell adhesion, migration and differentiation. This is one of the isozymes catalyzing that reaction (By similarity). [UniProt]
Calculated Mw	63 kDa
Cellular Localization	Membrane; Multi-pass membrane protein. [UniProt]

Images



ARG59584 anti-HAS3 antibody WB image

Western blot: 25 µg of Rat heart, Mouse heart and SKOV3 cell lysates stained with ARG59584 anti-HAS3 antibody at 1:1000 dilution.