

## ARG65829 anti-Follicular Dendritic Cells antibody [BU10] (FITC)

Package: 100 µl

Store at: 4°C

### Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [BU10] recognizes Follicular Dendritic Cells
Tested Reactivity	Hu
Tested Application	FACS
Specificity	Recognizes a molecule that is expressed on human follicular dendritic cells (FDCs).
Host	Mouse
Clonality	Monoclonal
Clone	BU10
Isotype	IgM
Target Name	Follicular Dendritic Cells
Species	Human
Immunogen	Follicular Dendritic Cells
Conjugation	FITC

### Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purified.
Buffer	PBS (pH 7.2), 0.09% Sodium azide and 1% BSA.
Preservative	0.09% Sodium azide
Stabilizer	1% BSA
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.