

## ARG65844 anti-Erythroid Cell / Ly76 antibody [TER-119] (PE)

Package: 100 µl  
Store at: 4°C

### Summary

Product Description	PE-conjugated Rat Monoclonal antibody [TER-119] recognizes Erythroid Cell / Ly76
Tested Reactivity	Ms
Tested Application	FACS
Specificity	The anti-mouse erythroid cell (Ly-76) mAb is selectively reactive with both fetal and adult erythroid cells. This monoclonal antibody (clone: TER-119) is specific for cells at stages from early proerythroblast to mature erythrocytes. TER-119 is reported to react with 20-25% of bone marrow cells and 2-3% of spleen cells but not with thymocytes or lymph node cells.
Host	Rat
Clonality	Monoclonal
Clone	TER-119
Isotype	IgG2b
Target Name	Erythroid Cell / Ly76
Species	Mouse
Immunogen	Mouse Erythroid Cell / Ly76
Conjugation	PE

### Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purified.
Buffer	PBS (pH 7.2), 0.09% Sodium azide and 1% BSA.
Preservative	0.09% Sodium azide
Stabilizer	1% BSA
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.