

ARG66303 Rabbit anti-Duck IgG (H+L) antibody (HRP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	HRP-conjugated Rabbit Polyclonal antibody recognizes Duck IgG (H+L)
Tested Reactivity	Duck
Tested Application	Dot, ELISA, ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Duck
Immunogen	Duck IgG
Conjugation	HRP

Application Instructions

Application table	Application	Dilution
	Dot	Assay-dependent
	ELISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified.
Buffer	PBS (pH 7.2)
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.