

ARG67206 anti-Feline Calicivirus VP1 Protein antibody [3C10]

Package: 100 µg
Store at: -20°C

Summary


Product Description	Mouse Monoclonal antibody [3C10] recognizes Feline Calicivirus VP1 Protein
Tested Reactivity	Virus
Tested Application	ICC/IF
Host	Mouse
Clonality	Monoclonal
Clone	3C10
Target Name	Feline Calicivirus VP1 Protein
Species	Virus
Immunogen	Full Virus of Feline Calicivirus
Conjugation	Un-conjugated
Alternate Names	Feline calicivirus, FCV

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:250 - 1:750
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS, 0.05% Sodium azide and 20% Glycerol.
Preservative	0.05% Sodium azide
Stabilizer	20% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG67206 anti-Feline Calicivirus VP1 Protein antibody [3C10] ICC/IF image

Immunofluorescence: ARG67206 anti-Feline Calicivirus VP1 Protein antibody [3C10] were used for detecting Feline Calicivirus VP1 Protein.