

ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198)
Tested Reactivity	Hu, Ms, Rat
Tested Application	ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	PKA C alpha + PKA C beta + PKA C gamma
Species	Human
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T198 of human PKA C alpha + PKA C beta + PKA C gamma.
Conjugation	Un-conjugated
Alternate Names	PRKACA; PKACA; cAMP-dependent protein kinase catalytic subunit alpha; PKA-alpha; PRKACB; cAMP-dependent protein kinase catalytic subunit beta; PKA-beta; PRKACG; cAMP-dependent protein kinase catalytic subunit gamma; PKA-gamma

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:100 - 1:500
	IHC-P	1:100 - 1:200
	WB	1:500 - 1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Observed Size	42 kDa	

Properties


Form	Liquid
Purification	Affinity purification with immunogen
Buffer	0.42% Potassium phosphate, 0.87% NaCl, 0.01% Sodium azide and 30% Glycerol.
Preservative	0.01% Sodium azide
Stabilizer	30% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot

and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note For laboratory research only, not for drug, diagnostic or other use.


Bioinformation

Gene Symbol	PRKACA; PRKACB; PRKACG
Gene Full Name	cAMP-dependent protein kinase catalytic subunit alpha / beta / gamma
Background	Protein kinase A (PKA, cAMP-dependent protein kinase) is a key element of a ubiquitous signaling pathway important in the cell cycle, cellular communication, memory formation and behavior. PKA is composed of two catalytic (PKAc; Protein Kinase A catalytic subunit) and two regulatory subunits (PKAr). Upon binding cAMP, the complex dissociates to PKAr dimer and two activated PKAc ser/thr protein kinase catalytic monomers. The released PKAc can translocate into the nucleus and exert a regulatory role in the activation of multiple nuclear hormone receptors. However, PKAc-mediated activation of tonicity-dependent gene expression is cAMP independent. Humans express three types of PKAc subunit – PKAc alpha is present in most human tissues, PKAc beta and gamma are tissue-specific, the later is found in testes.
Function	Phosphorylates a large number of substrates in the cytoplasm and the nucleus. Regulates the abundance of compartmentalized pools of its regulatory subunits through phosphorylation of PJA2 which binds and ubiquitinates these subunits, leading to their subsequent proteolysis. Phosphorylates CDC25B, ABL1, NFKB1, CLDN3, PSMC5/RPT6, PJA2, RYR2, RORA and VASP. RORA is activated by phosphorylation. Required for glucose-mediated adipogenic differentiation increase and osteogenic differentiation inhibition from osteoblasts. Involved in the regulation of platelets in response to thrombin and collagen; maintains circulating platelets in a resting state by phosphorylating proteins in numerous platelet inhibitory pathways when in complex with NF-kappa-B (NFKB1 and NFKB2) and I-kappa-B-alpha (NFKBIA), but thrombin and collagen disrupt these complexes and free active PRKACA stimulates platelets and leads to platelet aggregation by phosphorylating VASP. Prevents the antiproliferative and anti-invasive effects of alpha-difluoromethylornithine in breast cancer cells when activated. RYR2 channel activity is potentiated by phosphorylation in presence of luminal Ca(2+), leading to reduced amplitude and increased frequency of store overload-induced Ca(2+) release (SOICR) characterized by an increased rate of Ca(2+) release and propagation velocity of spontaneous Ca(2+) waves, despite reduced wave amplitude and resting cytosolic Ca(2+). PSMC5/RPT6 activation by phosphorylation stimulates proteasome. Negatively regulates tight junctions (TJs) in ovarian cancer cells via CLDN3 phosphorylation. NFKB1 phosphorylation promotes NF-kappa-B p50-p50 DNA binding. Involved in embryonic development by down-regulating the Hedgehog (Hh) signaling pathway that determines embryo pattern formation and morphogenesis. Prevents meiosis resumption in prophase-arrested oocytes via CDC25B inactivation by phosphorylation. May also regulate rapid eye movement (REM) sleep in the pedunculo pontine tegmental (PPT). Phosphorylates APOBEC3G and AICDA. Isoform 2 phosphorylates and activates ABL1 in sperm flagellum to promote spermatozoa capacitation. [UniProt]
Research Area	Cancer antibody; Metabolism antibody; Signaling Transduction antibody
Calculated Mw	40 kDa
PTM	Asn-3 is partially deaminated to Asp giving rise to 2 major isoelectric variants, called CB and CA respectively. Autophosphorylated. Phosphorylation is enhanced by vitamin K(2). Phosphorylated on threonine and serine residues. Phosphorylation on Thr-198 is required for full activity. Phosphorylated at Tyr-331 by activated receptor tyrosine kinases EGFR and PDGFR; this increases catalytic efficiency.




ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody IHC-P image

Immunohistochemistry: Human breast cancer stained with ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody.




ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody ICC/IF image

Immunofluorescence: HEK293T stained with ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody.




ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody WB image

Western blot: HeLa stained with ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody.



ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody WB image

Western blot: Rat testis stained with ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody.



ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody WB image

Western blot: Mouse testis stained with ARG67343 anti-PKA C alpha + PKA C beta + PKA C gamma phospho (Thr198) antibody.