

ARG67353 anti-MSK2 / RSK-B Phospho (Thr568) antibody

Package: 100 µl, 50 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes MSK2 / RSK-B Phospho (Thr568)
Tested Reactivity	Hu, Ms
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	MSK2 / RSK-B
Species	Human
Immunogen	The antiserum was produced against synthesized peptide derived from human MSK2 / RSK-B around the phosphorylation site of Thr568. AA range:531-580
Conjugation	Un-conjugated
Alternate Names	MSK2; Ribosomal protein S6 kinase alpha-4; S6K-alpha-4; 90 kDa ribosomal protein S6 kinase 4; Nuclear mitogen- and stress-activated protein kinase 2; Ribosomal protein kinase B; RSKB

Application Instructions

Application table	Application	Dilution
	IHC-P	1:100 - 1:300
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Observed Size	95 kDa	


Properties

Form	Liquid
Purification	Affinity purification with immunogen
Buffer	PBS, 0.02% Sodium azide, 50% Glycerol, and 0.5% BSA
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol, 0.5% BSA
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	RPS6KA4
Gene Full Name	Ribosomal protein S6 kinase alpha-4
Background	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates various substrates, including CREB1 and ATF1. The encoded protein can also phosphorylate histone H3 to regulate certain inflammatory genes. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]
Function	Serine/threonine-protein kinase that is required for the mitogen or stress-induced phosphorylation of the transcription factors CREB1 and ATF1 and for the regulation of the transcription factor RELA, and that contributes to gene activation by histone phosphorylation and functions in the regulation of inflammatory genes. Phosphorylates CREB1 and ATF1 in response to mitogenic or stress stimuli such as UV-C irradiation, epidermal growth factor (EGF) and anisomycin. Plays an essential role in the control of RELA transcriptional activity in response to TNF. Phosphorylates 'Ser-10' of histone H3 in response to mitogenics, stress stimuli and EGF, which results in the transcriptional activation of several immediate early genes, including proto-oncogenes c-fos/FOS and c-jun/JUN. May also phosphorylate 'Ser-28' of histone H3. Mediates the mitogen- and stress-induced phosphorylation of high mobility group protein 1 (HMGN1/HMG14). In lipopolysaccharide-stimulated primary macrophages, acts downstream of the Toll-like receptor TLR4 to limit the production of pro-inflammatory cytokines. Functions probably by inducing transcription of the MAP kinase phosphatase DUSP1 and the anti-inflammatory cytokine interleukin 10 (IL10), via CREB1 and ATF1 transcription factors. [UniProt]
Calculated Mw	84 kDa
PTM	Phosphoprotein. [UniProt]
Cellular Localization	Nucleus. [UniProt]

Images



ARG67353 anti-MSK2 / RSK-B Phospho (Thr568) antibody IHC-P image

Immunohistochemistry: Human brain, stained with ARG67353 anti-MSK2 / RSK-B Phospho (Thr568) antibody.



ARG67353 anti-MSK2 / RSK-B Phospho (Thr568) antibody WB image

Western blot: 293 stained with ARG67353 anti-MSK2 / RSK-B Phospho (Thr568) antibody.