

ARG67354 anti-PTH LH antibody

Package: 100 µl, 50 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes PTH LH
Tested Reactivity	Hu, Ms
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	PTH LH
Immunogen	Synthesized peptide derived from the Internal region of human PTH LH. AA range: 101-150
Conjugation	Un-conjugated
Alternate Names	Parathyroid hormone-related protein; PTH-rP; PLP; PTHrP[107-139]; PTHRP; Parathyroid hormone-like protein; PTHR; PTHrP; HHM; BDE2

Application Instructions

Application table	Application	Dilution
	IHC-P	1:50 - 1:300
	WB	1:500 -1:2000
Application Note	* For IHC, we suggest antigen retrieval with TE buffer pH 9.0 * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Affinity purification with immunogen
Buffer	PBS, 0.02% Sodium azide, 50% Glycerol, and 0.5% BSA
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol, 0.5% BSA
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation


Gene Symbol	PTHLH
Gene Full Name	parathyroid hormone-like hormone
Background	The protein encoded by this gene is a member of the parathyroid hormone family. This hormone, via its receptor, PTHR1, regulates endochondral bone development and epithelial-mesenchymal interactions during the formation of the mammary glands and teeth. It is responsible for most cases of humoral hypercalcemia of malignancy, and mutations in this gene are associated with brachydactyly type E2 (BDE2). Alternatively spliced transcript variants have been found for this gene. There is also evidence for alternative translation initiation from non-AUG (CUG and GUG) start sites, downstream of the initiator AUG codon, resulting in nuclear forms of this hormone. [provided by RefSeq, Nov 2013]
Function	Neuroendocrine peptide which is a critical regulator of cellular and organ growth, development, migration, differentiation and survival and of epithelial calcium ion transport. Regulates endochondral bone development and epithelial-mesenchymal interactions during the formation of the mammary glands and teeth. Required for skeletal homeostasis. Promotes mammary mesenchyme differentiation and bud outgrowth by modulating mesenchymal cell responsiveness to BMPs. Upregulates BMPR1A expression in the mammary mesenchyme and this increases the sensitivity of these cells to BMPs and allows them to respond to BMP4 in a paracrine and/or autocrine fashion. BMP4 signaling in the mesenchyme, in turn, triggers epithelial outgrowth and augments MSX2 expression, which causes the mammary mesenchyme to inhibit hair follicle formation within the nipple sheath (By similarity). Promotes colon cancer cell migration and invasion in an integrin alpha-6/beta-1-dependent manner through activation of Rac1. Osteostatin is a potent inhibitor of osteoclastic bone resorption. [UniProt]
Calculated Mw	20 kDa
PTM	There are 3 principal secretory forms, called PTHrP[1-36], PTHrP[38-94], and osteostatin (PTHrP[107-139]) arising from endoproteolytic cleavage of the initial translation product. Each of these secretory forms is believed to have one or more of its own receptors that mediates the normal paracrine, autocrine and endocrine actions.

Images



ARG67354 anti-PTHLH antibody IHC-P image

Immunohistochemistry: Human cervical carcinoma stained with ARG67354 anti-PTHLH antibody at 1:200 dilution.



ARG67354 anti-PTHLH antibody WB image

Western blot: HeLa stained with ARG67354 anti-PTHLH antibody at 1:800 dilution.