


ARG70077 Human IL33 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL33 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL33
Species	Human
A.A. Sequence	Ser112 - Thr270
Expression System	E. coli
Activity	Active
Activity Note	

Images



ARG70077 Human IL33 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70077 Human IL33 recombinant protein (Active) (His-tagged, C-ter).