

## ARG70131 Human IGF2 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human IGF2 recombinant protein
Tested Application	SDS-PAGE
Target Name	IGF2
Species	Human
A.A. Sequence	Ala25 - Glu91
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce MCF-7 cells proliferation. The ED50 for this effect is $3 \times 10^5$ IU/mg.
Alternate Names	Somatomedin-A; Insulin-like growth factor II; IGF-II; PP9974; GRDF; C11orf43; T3M-11-derived growth factor

### Properties

Form	Powder
Purification Note	Endotoxin level is 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

### Images



ARG70131 Human IGF2 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70131 Human IGF2 recombinant protein (Active) (His-tagged, N-ter).