

ARG70217 SARS-CoV-2 nucleocapsid recombinant protein (Tag Free)

Package: 100 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, SARS-CoV-2 nucleocapsid recombinant protein
Tested Application	SDS-PAGE
Target Name	SARS-CoV-2 nucleocapsid protein
Species	Virus
A.A. Sequence	Met1 - Ala419
Expression System	E. coli

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS containing 10 mM potassium carbonate (pH 7.4)
Reconstitution	

Images



SARS-CoV-2 nucleocapsid
recombinant protein

ARG70217 SARS-CoV-2 nucleocapsid recombinant protein (Tag Free)
SDS-PAGE image

SDS-PAGE analysis of ARG70217 SARS-CoV-2 nucleocapsid
recombinant protein (Tag Free).