

ARG70766 Mouse CD72 recombinant protein (His-tagged, C-ter)

Package: 20 µg
Store at: -20°C

Summary


Product Description	CHO expressed, His-tagged (C-ter) CD72 recombinant protein
Tested Application	SDS-PAGE
Target Name	CD72
Species	Mouse
A.A. Sequence	Arg117 - Asn354
Expression System	CHO
Alternate Names	CD antigen CD72; LYB2; B-cell differentiation antigen CD72; Lyb-2; CD72b

Properties

Form	Powder
Purification Note	Endotoxin level is less than 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 90% (by SDS-PAGE); > 95% (by HPLC)
Buffer	PBS (pH 7.4)
Reconstitution	It is recommended to reconstitute the lyophilized protein in sterile water to a concentration not less than 200 µg/mL and incubate the stock solution for at least 20 min at room temperature to make sure the protein is dissolved completely.
Storage instruction	For long term, lyophilized protein should be stored at -20°C or -80°C. After reconstitution, aliquot and store at -20°C or -80°C for up to one month. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	CD72
Gene Full Name	CD72 antigen
Function	Plays a role in B-cell proliferation and differentiation. [UniProt]
PTM	Phosphorylated upon engagement of the B-cell receptor, probably by LYN or SYK. Phosphorylation at Tyr-7 is important for interaction with PTPN6/SHP-1 (By similarity).



ARG70766 Mouse CD72 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70766 Mouse CD72 recombinant protein (His-tagged, C-ter).