

## ARG80369 Human Anti-Spermatozoa antibody (ASA) ELISA Kit

Package: 96 wells

Store at: 4°C

### Summary

Product Description	<p>ARG80369 Human Anti-Spermatozoa antibody (ASA) ELISA Kit is a quick, reliable, quantitative test for the detection of antibodies directed against human spermatozoa. This test may be used with serum.</p> <p>Please note: the terms anti-spermatozoa antibodies , anti-sperm antibodies and sperm antibodies are equivalent. In these instructions the rather unwieldy but correct term anti-spermatozoa antibodies is used.</p>
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Anti-Spermatozoa antibody (ASA)
Sensitivity	25 U/ml
Sample Type	Serum.
Standard Range	31 - 250 U/ml
Sample Volume	50 µl

### Application Instructions

Assay Time	1 h, 1 h (37°C), 30 min (RT)
------------	------------------------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Resrarch Area	Immune System kit
---------------	-------------------