

ARG80384 Human dsDNA Antibody Screen ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80384 Human dsDNA Antibody Screen ELISA Kit is an Enzyme Immunoassay kit for the quantification of IgG, IgM and IgA autoantibodies to dsDNA in human serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	ds-DNA Screen
Sensitivity	1 U/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	12.5 - 200 U/ml
Sample Volume	10 µl

Application Instructions

Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: DNA antibodies ; DNA ELISA Kits ; DNA Duos / Panels ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Immune System kit