

ARG80403 Human Phosphatidyl Inositol antibody IgG / IgM ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80403 Human Phosphatidyl Inositol antibody IgG / IgM ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Phosphatidyl Inositol Ab IgG / IgM in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Phosphatidyl Inositol antibody IgG / IgM
Sensitivity	0.5 U/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	6.3 - 100 U/ml
Sample Volume	10 µl

Application Instructions

Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Cell Biology and Cellular Response kit; Immune System kit; Metabolism kit