

ARG80403 Human Phosphatidyl Inositol antibody IgG / IgM ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

| | |
|---------------------|---|
| Product Description | ARG80403 Human Phosphatidyl Inositol antibody IgG / IgM ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Phosphatidyl Inositol Ab IgG / IgM in serum and plasma (EDTA, heparin, citrate). |
| Tested Reactivity | Hu |
| Tested Application | ELISA |
| Target Name | Phosphatidyl Inositol antibody IgG / IgM |
| Sensitivity | 0.5 U/ml |
| Sample Type | Serum and plasma (EDTA, heparin, citrate). |
| Standard Range | 6.3 - 100 U/ml |
| Sample Volume | 10 µl |

Application Instructions

| | |
|------------|------------------------|
| Assay Time | 30, 15, 15, 5 min (RT) |
|------------|------------------------|

Properties

| | |
|---------------------|--|
| Form | 96 well |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |

Bioinformation

| | |
|---------------|---|
| Resrarch Area | Cell Biology and Cellular Response kit; Immune System kit; Metabolism kit |
|---------------|---|