

ARG80405 Human Phospholipid Screen IgG / IgM ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG80405 Human Phospholipid Screen IgG / IgM ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Phospholipid Screen IgG / IgM in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Phospholipid Screen IgG / IgM
Sensitivity	0.5 U/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	6.3 - 100 U/ml
Sample Volume	10 µl

Application Instructions

Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Resrarch Area	Cell Biology and Cellular Response kit; Immune System kit; Metabolism kit
---------------	---