

ARG80416 Human Scl70 antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80416 Human Scl70 antibody ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Scl70 Ab in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Scl70 antibody
Conjugation	HRP
Sensitivity	1 U/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	12.5 - 200 U/ml
Sample Volume	10 µl

Application Instructions


Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Scl70 antibody ELISA Kits: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Immune System kit



ARG80416 Human Scl70 antibody ELISA Kit standard curve image

ARG80416 Human Scl70 antibody ELISA Kit results of a typical standard run with optical density reading at 450 / 620 nm.