

ARG80467 Human Metanephrine ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80467 Human Metanephrine ELISA Kit is an Enzyme Immunoassay kit for the quantification of Metanephrine in urine.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Metanephrine
Sensitivity	13 ng/ml
Sample Type	Urine.
Standard Range	20 - 2000 ng/ml
Sample Volume	25 µl

Application Instructions


Assay Time	30, 15, 15 min (RT/shaker)
------------	----------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Resrarch Area	Metabolism kit; Neuroscience kit; Signaling Transduction kit
---------------	--



ARG80467 Human Metanephrine ELISA Kit standard curve example image

The representative standard curve of ARG80467 Human Metanephrine ELISA Kit. The standard curve is for demonstration only and cannot be used in place of data generations at the time of assay. The standard curve should be generated each time the assay is performed.