

ARG80588 Human Mycobacterium tuberculosis IgG antibody ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG80588 Human Mycobacterium tuberculosis IgG antibody ELISA Kit has been designed for the determination of specific IgG antibodies against Mycobacterium tuberculosis in serum and plasma (EDTA, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Mycobacterium tuberculosis IgG antibody
Sensitivity	1.09 U/ml
Sample Type	Serum and plasma (EDTA, heparin).
Standard Range	1 - 150 U/ml
Sample Volume	5 µl

Application Instructions


Assay Time	1 h, 30, 20 min
------------	-----------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.


Bioinformation

Resrarch Area	Microbiology and Infectious Disease kit
---------------	---



ARG80588 Human Mycobacterium tuberculosis IgG antibody ELISA Kit standard curve image

ARG80588 Human Mycobacterium tuberculosis IgG antibody ELISA Kit results of a typical standard.



ARG80588 Human Mycobacterium tuberculosis IgG antibody ELISA Kit standard curve image

ARG80588 Human Mycobacterium tuberculosis IgG antibody ELISA Kit Results of a typical standard run with optical density reading at 450 nm.