

## ARG80663 Canine TSH ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG80663 Canine TSH ELISA Kit is an enzyme immunoassay kit for the quantification of canine TSH (thyrotropin) in serum and plasma (EDTA).
Tested Reactivity	Dog
Tested Application	ELISA
Target Name	TSH
Conjugation	HRP
Sensitivity	0.01 ng/ml
Sample Type	Serum and plasma (EDTA).
Standard Range	0.2 - 5.2 ng/ml
Sample Volume	100 µl
Alternate Names	TSH-B; TSH-beta; Thyroid-stimulating hormone subunit beta; Thyrotropin subunit beta; Thyrotropin beta chain; Thyrotropin alfa; TSH-BETA

### Application Instructions

Assay Time	2 h (shaking), 30 min
------------	-----------------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Gene Symbol	TSHB
Gene Full Name	thyroid stimulating hormone, beta
Background	The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological specificity. Thyroid stimulating hormone functions in the control of thyroid structure and metabolism. The protein encoded by this gene is the beta subunit of thyroid stimulating hormone. Mutations in this gene are associated with congenital central and secondary hypothyroidism and Hashimoto's thyroiditis. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, May 2013]
Function	Indispensable for the control of thyroid structure and metabolism. [UniProt]


Highlight

Related products:  
[TSH antibodies](#); [TSH ELISA Kits](#); [TSH Duos / Panels](#);  
New ELISA data calculation tool:  
[Simplify the ELISA analysis by GainData](#)

Research Area

Cancer kit; Metabolism kit; Neuroscience kit; Signaling Transduction kit

Images



ARG80663 Canine TSH ELISA Kit standard curve image  
ARG80663 Canine TSH ELISA Kit results of a typical standard.